



Germo[®]
S.P.A.

Sede Legale Stabilimento e Uffici: 20032 CORMANO (Milano) Via Giotto, 19/21
Telefono +39-0266301938 r.a. - Telefax +39-0266301939
E-mail: info@germodis.com www.germo.it

Capitale Sociale € 156.000,00 interamente versato - R.E.A. 497966 Cod. Fisc. e Part. IVA 00772350153 - Iscr. 125750 R.I. Canc. Trib. MI

ANTIOXIDANT LINE - TECHNICAL DATASHEET

ANTIOSSIDANTE

Sodium Nitrite in sachets with antioxidant action.

Antioxidant with Sodium nitrite to be added to disinfectant or sterilizer solutions to prevent the oxidation of surgical instruments. Especially indicated in case of long contact times.



INSTRUCTIONS FOR USE

Dissolve a sachet of Sodium Nitrite Antioxidant in the solution with the disinfectant / sterilizer and follow the instructions on the label of the product used.

PACKAGING

Internal Code:	Primary packaging:	Secondary packaging:
Q720	Single sachet	Box containing 500 sachets



Germo[®]
S.P.A.

Sede Legale Stabilimento e Uffici: 20032 CORMANO (Milano) Via Giotto, 19/21
Telefono +39-0266301938 r.a. - Telefax +39-0266301939
E-mail: info@germodis.com www.germo.it

Capitale Sociale € 156.000,00 interamente versato - R.E.A. 497966 Cod. Fisc. e Part. IVA 00772350153 - Iscr. 125750 R.I. Canc. Trib. MI

ANTIOXIDANT LINE - TECHNICAL DATASHEET

ANTIOSSIDANTE

COMPOSITION

Sodium Nitrite at 99%.

PRODUCT FORM

White to slightly pale crystalline powder.

HOW TO STORAGE

Store in a cool, dry place, away from heat sources.

QUALITY CONTROL

The components (raw materials, containers, labels, etc.) and the processing steps of each lot are regularly and carefully checked internally following the procedures of the Organization Quality Management System certified UNI EN ISO 9001.

CONFIDENTIAL INFORMATION FOR PROFESSIONAL USERS